



positioning
patients for life.®

SofTouch™

Knee Wedges & Foot Block

SofTouch™ Knee Wedges & Foot Block

The SofTouch™ Knee Wedges & Foot Block are designed to provide repeatable positioning and are contoured to relieve pressure point stress and improve patient comfort. They are lightweight and easy to clean.

Features:

- Ergonomically designed for patient comfort while maintaining reproducible positioning
- Indexes to locating bar in 35 mm increments
- Centering groove to aid in laser alignment
- Optional Clip-on leg and foot straps for added immobilization
- Cut-out for ERB syringe
- Lightweight
- Easy to clean



Qfix
Avondale, PA
+1 610.268.0585
sales@Qfix.com
www.Qfix.com

RT-4489KW-01	SofTouch Knee Wedge - Small (includes strap)
RT-4489KW	SofTouch Knee Wedge - Large (includes strap)
RT-4489FB	SofTouch Foot Blocks (includes strap)