

Stradivarius™

Compression Belt

The Stradivarius™ Compression Belt is an easy to use abdominal compression system for managing tumor motion. The belt is fastened around the patient below the rib cage, so that the pressure plate is applied to the diaphragm. Using the compact, one-body design manometer and inflation bulb, the bladder is inflated, suppressing the diaphragm and reducing diaphragm excursions. Two different size belts make it highly adaptable to patients of all shapes and sizes. The Stradivarius™ Compression Belt is used in conjunction with the Stradivarius™ SBRT system and vacuum cushion that is custom fitted for complete body immobilization. It can also be used as a standalone device.

Product Highlights:

- Easy to use abdominal motion management system allows for compression and suppressed tumor motion
- Compact, one-body design manometer has an operational range of 0-300 mmHg and is designed to monitor delivered pressures
- Two belt sizes (115 cm / 150 cm) accommodate every patient of every size



Order Information:

RT-4558CB01 | Stradivarius™ Compression Belt