



positioning
patients for life.®

**Lightweight
Only 2.5 kg!**

Portrait™


Intracranial, Head & Neck Device

Stable, reliable positioning

The lightweight Portrait™ Intracranial, Head and Neck Device provides superior immobilization for a variety of brain, head and neck radiotherapy treatments.




Portrait Intracranial,
Head & Neck Device
RT-4552



Portrait MR Intracranial,
Head & Neck Device
RT-4552MR

Features of Portrait Head & Neck Device

- Easy to use
- Lightweight
- Daily reproducibility



kVue Portrait Insert
RT-4552KV

The Portrait™ Intracranial, Head & Neck Device and Portrait S-Frame Thermoplastic Masks enable radiotherapy treatment with confidence, providing a high level of stability by limiting patient motion and reducing the probability of positioning errors.

The device indexes to the couch top with a two pin Locating Bar or can utilize Direct Indexing™ for increased reproducibility and stability. The Portrait accommodates a variety of thermoplastic masks and head supports* and accommodates supine and prone setups for patient specific treatments. The Portrait is rigid while providing minimal attenuation making it ideal for full 360-degree IMRT treatments.

The Portrait™ MR Intracranial, Head and Neck Device is CT and MR compatible for enhanced imaging and can also be used for treatment.



RT-4552	Portrait Head & Neck Device
RT-4552MR	Portrait MR Head & Neck Device
RT-4552KV	kVue Portrait Insert

* Requires a Silverman Head Support. Additional accessories are also available.



Qfix
Avondale, PA
+1 610.268.0585
sales@Qfix.com
www.Qfix.com