

# Reproducible Positioning & Patient Comfort

NEW



**FLEXLOCK™**

**Knee & Foot Supports**



The sophisticated functionality of the FlexLock™ Knee and Foot Supports provide customizable indexing options for positioning the lower abdomen, pelvis, and lower extremities. These versatile supports index directly to MRI, CT, and treatment couch tops with standard two-pin locating bars and allow natural, comfortable positioning of the legs and feet.

The novel FlexLock™ Knee Support features integrated height adjustability so that only one device is required for a variety of patients. Superior to inferior adjustments can be made while the patient is positioned to fine tune the indexed location and improve patient comfort. This can reduce the need to remove the FlexLock™ Knee Support just to reposition the locating bar and can expedite the process of making minor adjustments.

The FlexLock™ Foot Support can be angled to accommodate patient comfort and offers superior to inferior adjustments with the use of the locating bar on the imaging or treatment surface to maintain a reproducible setup position.

## Features:

- Ergonomically designed for patient comfort while maintaining reproducible positioning
- Optional Clip-on leg and foot straps for added immobilization
- Highly durable and easily cleanable

<b>RT-4800-FS</b>	FlexLock™ Foot Support
<b>RT-4800-KS</b>	FlexLock™ Knee Support
<b>RT-4489STRAP</b>	Leg & Foot Strap Assembly, Set of 2



Qfix  
Avondale, PA  
+1 610.268.0585  
sales@Qfix.com  
www.Qfix.com

2008636\_EN\_B

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides. ©2020 Qfix. All rights reserved. Qfix and FlexLock are trademarks of Qfix.