



positioning
patients for life.®

TECHNICAL DATA SHEET

RT-4546

Access™ Supine Breast & Lung




Product Name	Access Supine Breast & Lung
Product Number	RT-4546
General Description	Patient Positioning System
Indications for Use	This device is intended to immobilize, position and reposition patients undergoing radiation therapy

TECHNICAL DATA SHEET – 2009049_A

Technical Details

Materials	Foam Cored Composite
MR Compatibility	NOT MR Compatible
Weight (Approx.)	13 lbs (5.9 kg) - with no accessories
Patient Weight Capacity	550 lbs (249 kg)
Dosimetric Properties	Water Equivalence: ~ 3 mm @ 6 MV in the grid pattern regions The structural regions of the Access Supine have a water equivalence of 5 mm at 6MV

Technical Drawing



Access™ Supine Breast & Lung
Access is a trademark of Qfix.

TECHNICAL DATA SHEET – 2009049_A