



positioning
patients for life.®

TECHNICAL DATA SHEET

RT-4525L

AccuFix™ Head & Neck Device




Product Name	AccuFix Head & Neck Device
Product Number	RT-4525L
General Description	Patient Positioning System
Indications for Use	This device is intended to immobilize, position and reposition patients undergoing radiation therapy

TECHNICAL DATA SHEET – 2009050_A

Technical Details

Materials	Foam Cored Composite
MR Compatibility	NOT MR Compatible
Weight (Approx.)	5.5 kg (12 lbs)
Patient Weight Capacity	249 kg (550 lbs)
Dosimetric Properties	Aluminum Equivalence: 0.6 mm @ 100 kVp

Technical Drawing



AccuFix™ Head & Neck Device
AccuFix is a trademark of Qfix

TECHNICAL DATA SHEET – 2009050_A