



positioning  
patients for life.®

# TECHNICAL DATA SHEET

## RT-4800-01

### Alta™ Multipurpose Device



<b>Product Name</b>	Alta Multipurpose Device
<b>Product Number</b>	RT-4800-01
<b>General Description</b>	Patient Positioning System
<b>Indications for Use</b>	The Alta™ Multipurpose Device is intended to immobilize, support, position, and transfer patients undergoing radiotherapy procedures using electrons, photons, protons, including SBRT and SRS; imaging procedures such as X-ray, computed tomography, and magnetic resonance imaging; and other procedures involving transfer of a patient.


**TECHNICAL DATA SHEET – 2008432\_C**

The Alta™ Multipurpose Device is designed to interface with the Symphony™ Patient Transport System and other positioning devices, such as couch tops, thermoplastic masks, setup-specific treatment devices and adapters, positioning supports and pads.

# Technical Details

<b>Materials</b>	Foam Cored Composite
<b>MR Compatibility</b>	MR Safe
<b>Weight (Approx.)</b>	5.4 kg (12 lbs)
<b>Patient Weight Capacity</b>	226 kg (500 lbs)
<b>Dosimetric Properties</b>	At <b>green</b> marker: Water Equivalence: 4.7 mm at 6 MV (10 cm x 10 cm field size) Aluminum Equivalence: 0.73 mm at 100 kVp

## Technical Drawing



Alta™ Multipurpose Device  
Alta is a trademark of Qfix.

### TECHNICAL DATA SHEET – 2008432\_C