




positioning
patients for life.®

TECHNICAL DATA SHEET

RT-4588-01

DoseMax™ Standalone Device




Product Name	DoseMax Standalone Device
Product Number	RT-4588-01
General Description	Patient Positioning System
Indications for Use	This device is intended to immobilize, position and reposition patients undergoing radiation therapy

TECHNICAL DATA SHEET – 2009052_A

Technical Details

Materials	Cored Composite	
MR Compatibility	NOT MR Compatible	
Weight (Approx.)	6.5 kg (14 lbs)	
Patient Weight Capacity	249 kg (550 lbs)	
Dosimetric Properties	Water Equivalence:	2.0 mm @ 6 MV

Technical Drawing



DoseMax™ Standalone Device
DoseMax is a trademark of Qfix.

TECHNICAL DATA SHEET – 2009052_A