




positioning
patients for life.®

TECHNICAL DATA SHEET

RT-4551KV18

kVue™ DoseMax™ Insert




Product Name	kVue DoseMax Insert
Product Number	RT-4551KV18
General Description	Patient Positioning System
Indications for Use	This device is intended to immobilize, position and reposition patients undergoing radiation therapy

TECHNICAL DATA SHEET – 2009058_A

Technical Details

Materials	Cored Composite	
MR Compatibility	NOT MR Compatible	
Weight (Approx.)	4 kg (9 lbs)	
Patient Weight Capacity	249 kg (550 lbs)	
Dosimetric Properties	Water Equivalence:	2.0 mm @ 6 MV

Technical Drawing



kVue™ DoseMax™ Insert
DoseMax and kVue are trademarks of Qfix

TECHNICAL DATA SHEET – 2009058_A