


One Simple Solution for a Variety of Treatment Setups

NEW kVue™ One Couch Top



RT-4551KV-01

*Shown with RT-4551KV19
kVue™ Low-e™ Insert*




Without moving rails

The kVue™ One Couch Top is a one-piece modular design with a homogenous composite base and no moving rails. It offers the ability to interchange treatment specific kVue™ inserts that mount directly into the top surface of the couch through the One Touch™ Latch. The light weight kVue™ One is compatible with versions of 6 Degrees of Freedom (DOF) linear accelerators and robotic bases. (5 mm WET without kVue™ Insert)

NEW kVue™ Low-e™ Insert

Low_{e-}

CALYPSO  OK




RT-4551KV19

*Shown with kVue™ SBRT Bridges
RT-4558-B1 (small),
RT-4558-B2 (large)*



The kVue™ Low-e™ Insert is a low electron generating, homogenous insert for enhanced imaging and treatment. It is the optimal choice for a wide range of radiotherapy treatments. (3.8 mm WET)

SBRT Accessories



Shown with low profile compression paddle

kVue™ SBRT Bridges

RT-4558-B1 (small)

RT-4558-B2 (large)

- Low profile bridges fit into most CT bore sizes
- Horizontal scale to aid in laser alignment
- Two sizes, Small (min: 30 cm & max: 38 cm)* and Large (min: 39 cm & max: 48 cm)* are completely adjustable in height

*Measurements are from the treatment surface to the underside of the bridge.



SofTouch™ Knee Wedges (large & small) & Foot Block

RT-4489KW (large)

RT-4489KW-01 (small)

& RT-4489FB

- Ergonomic design
- Indexes to locating bar in 35 mm increments
- Clip-on leg straps for added immobilization (included)
- Centering groove to aid in laser alignment



Shown with standard compression paddle

Compression Paddle

RT-4558KV4 (standard)

RT-4558KV4-01 (low profile)


- Quick release button
- Tactile compression
- Incremental lateral placement along bridge



ArmShuttle™ Elite

RT-4546C

- Lightweight
- Ergonomic design
- MR Safe & CT compatible




2-Chamber VacQfix™ Cushion

RT-4517-0002

3-Chamber VacQfix™ Cushion

RT-4517-0003


- Can be formed by one person
- Customizable full body immobilization
- Flared top for arm immobilization



VacQfix™ Indexer

RT-4517-IND01

- Quick, easy vacuum cushion indexing
- Adaptable for a variety of setups



Thigh Compression

VacQfix™ Cushion

RT-4517-4079VF15

- Complete leg immobilization
- Adjustable along bridge
- Snap-on attachment



ZiFix™, Abdominal Thoracic Motion Control System

RT-4558CB02

- Suppresses motion of lungs and liver for imaging & radiotherapy
- Compact, one-body design
- Three belt sizes

2008278_EN_A

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides.



Qfix
Avondale, PA
+1 610.268.0585
sales@Qfix.com
www.Qfix.com